

Formalin Dispensing Practices in Hospital Theatre and Pathology Settings.

A Technical White Paper for Facility Managers & Safety Leads

AUTHOR: Rebecca Burnett, Cabinet Technical Specialist

CONTENTS

- EXECUTIVE SUMMARY
- 1. FORMALIN USE IN HOSPITAL THEATRE AND PATHOLOGY SETTINGS
 - 1.1 Typical Use Cases
 - 1.2 Occupational Exposure Considerations
- 2. MANUAL FORMALIN DECANTING PROCESS CHARACTERISTICS
 - 2.1 Open Handling and Vapour Release
 - 2.2 Operator Dependency and Process Variability
 - 2.3 Position Within the Hierarchy of Control
- 3. CONTROLLED FORMALIN DISPENSING PROCESS CHARACTERISTICS
 - 3.1 Enclosed Dispensing Environments
 - 3.2 Air Management and Vapour Capture
 - 3.3 Process Consistency and Repeatability
- 4. EXPOSURE PATHWAY CONSIDERATIONS ACROSS DISPENSING METHODS
- 5. ALIGNMENT WITH THE HIERARCHY OF CONTROL
-
- 6. IMPLICATIONS FOR COMPLIANCE, COSHH ASSESSMENT, AND AUDIT
 - 6.1 Demonstrating Reasonably Practicable Controls
 - 6.2 Audit Readiness and Documentation
- 7. INTEGRATION INTO THEATRE AND PATHOLOGY WORKFLOWS
- 8. HUMAN FACTORS AND STAFF CONFIDENCE
- 9. EARLY ADOPTION AND FUTURE CONSIDERATIONS
- 10. CONCLUSION
- 11. GLOSSARY
- 12. STANDARDS DISCLAIMER
- APPENDIX A: EXAMPLE OF A CONTROLLED FORMALIN DISPENSING IMPLEMENTATION

EXECUTIVE SUMMARY

Formalin is a fundamental chemical in theatre and pathology workflows, widely used for tissue fixation and specimen preservation. Despite its utility, formalin contains formaldehyde, a recognised hazardous substance with known vapour exposure risks.

As healthcare environments place growing emphasis on exposure reduction, process consistency, and demonstrable compliance, controlled dispensing methods are being considered as part of wider risk mitigation strategies.

This white paper examines formalin dispensing practices from a process and risk-control perspective, outlining the attributes and implications of manual and controlled dispensing approaches, and offering insight into how controlled systems might support organised exposure management and compliance.



1. FORMALIN USE IN HOSPITAL THEATRE AND PATHOLOGY SETTINGS

1.1 Typical Use Cases

Formalin is commonly employed in:

- Theatre support workflows immediately following surgical procedures
- Pathology and histology departments for tissue fixation
- Back-bench specimen preparation prior to laboratory analysis

These activities often involve multiple transfers of formalin, across various work areas and by multiple staff members, reinforcing the importance of robust process controls.

1.2 Occupational Exposure Considerations

The principal health concern associated with formalin handling is inhalation of formaldehyde vapour. Exposure can occur during container opening, transfer, and spill events.

Cumulative exposure, even at low levels, can contribute to discomfort and potential long-term health effects, making exposure control an important element in occupational health strategies.

2. MANUAL FORMALIN DECANTING PROCESS CHARACTERISTICS

2.1 Open Handling and Vapour Release

Manual decanting involves transferring formalin from a bulk container into secondary vessels by hand. This process requires opening primary containers and pouring liquid, which inherently exposes operators to formaldehyde vapour. Each stage [lid removal, pour, reseal] represents a potential vapour release point.

2.2 Operator Dependency and Process Variability

The effectiveness of manual decanting in controlling vapour exposure depends on individual technique. Variability in pouring height, speed, and alignment can influence the degree of vapour generation and spill risk. This operator dependency makes it difficult to achieve consistent control outcomes across staff and shifts.

2.3 Position within the Hierarchy of Control

In the recognised Hierarchy of Control, manual decanting relies primarily on administrative controls (procedures, training) and PPE. While important, these measures do not address vapour generation at source and are inherently susceptible to human factors.

3. CONTROLLED FORMALIN DISPENSING PROCESS CHARACTERISTICS

3.1 Enclosed Dispensing Environments


Controlled dispensing systems are designed to contain the formalin transfer process within a defined, enclosed workspace. By limiting open handling, these systems reduce the potential for formaldehyde vapour to escape into the surrounding environment.

3.2 Air Management and Vapour Capture

Many controlled dispensing solutions incorporate vapour management features, such as activated carbon filtration or localised airflow control, which capture and remove airborne contaminants during dispensing. This approach helps manage ambient vapour concentrations and supports a safer handling environment.

3.3 Process Consistency and Repeatability

By standardising mechanical dispensing steps, controlled systems reduce process variability between operators. Measured dispensing and repeatable transfer mechanisms support consistent outcomes, which can strengthen procedural predictability and reduce reliance on individual technique.



4. EXPOSURE PATHWAY CONSIDERATIONS ACROSS DISPENSING METHODS

Exposure pathways differ between manual and controlled dispensing methods. Manual open pouring creates multiple uncontrolled release points, whereas enclosed dispensing focuses transfer within a managed environment.

Reducing exposure opportunities in this way aligns with COSHH expectations to control risk so far as is reasonably practicable.

5. ALIGNMENT WITH THE HIERARCHY OF CONTROL

The Hierarchy of Control prioritises measures that eliminate or reduce risk at source over those that rely on behavioural compliance.

Controlled dispensing represents an engineering control, directly addressing potential exposure pathways, whereas manual decanting remains dependent on lower-order procedural and PPE controls.

6. IMPLICATIONS FOR COMPLIANCE, COSHH ASSESSMENT, AND AUDIT

6.1 Demonstrating Reasonably Practicable Controls

COSHH assessments require employers to demonstrate that exposure to hazardous substances has been reduced as far as reasonably practicable. Incorporating engineering controls into formalin handling processes strengthens the evidence base for such assessments.

6.2 Audit Readiness and Documentation

Standardised, controlled dispensing processes support clearer documentation and training records, simplifying internal and external audits.

7. INTEGRATION INTO THEATRE AND PATHOLOGY WORKFLOWS

Controlled dispensing systems can be integrated into existing workflows with minimal disruption when selected appropriately.

Considerations include space constraints, compatibility with existing containers, maintenance requirements, and staff training.

8. HUMAN FACTORS AND STAFF CONFIDENCE

Controlled systems reduce reliance on manual pouring and operator technique, supporting safer routine practice and improving staff confidence when handling hazardous substances.

9. EARLY ADOPTION AND FUTURE CONSIDERATIONS

As exposure reduction expectations continue to evolve, early review of formalin handling practices can support long-term compliance and safer working environments. Controlled dispensing approaches may form part of this progression.

10. CONCLUSION

Evaluation of formalin dispensing practices demonstrates that transfer methodology materially influences exposure pathways, process variability, and the overall effectiveness of risk control measures.

Manual decanting introduces multiple potential vapour release points and relies predominantly on administrative controls and PPE, both of which are inherently operator-dependent.

Controlled dispensing approaches, by contrast, apply engineering principles to reduce exposure at source and support greater process consistency.

This alignment with higher-order measures within the Hierarchy of Control strengthens the evidential basis for COSHH assessments and compliance documentation.

Reviewing formalin dispensing methodology should therefore be considered a component of broader occupational risk management within theatre and pathology environments.

Selection of appropriate control measures must be informed by site-specific assessment, workflow requirements, and regulatory obligations, with the objective of achieving demonstrable and reasonably practicable exposure reduction.



11. GLOSSARY

- **COSHH**

Control of Substances Hazardous to Health Regulations

- **Engineering Control**

A control measure that physically reduces exposure

- **Exposure Pathway**

A route by which a hazardous substance may contact a person

- **Formalin**

An aqueous solution of formaldehyde used for tissue fixation

- **Hierarchy of Control**

Framework ranking risk controls by effectiveness

12. STANDARDS DISCLAIMER

This document references general UK health and safety principles, including COSHH and HSE guidance.

It is intended for informational purposes only and does not replace site-specific risk assessments or regulatory consultation.

APPENDIX A: EXAMPLE OF A CONTROLLED FORMALIN DISPENSING IMPLEMENTATION

Controlled dispensing systems are commercially available and may be specified to support exposure reduction objectives in hospital theatre and pathology environments.

These systems are designed to manage formalin transfer within an enclosed workspace, reducing reliance on manual decanting and open handling.

One example of such an implementation is the **Monmouth Scientific Circulaire® Chemical Dispensing Station**, an enclosed dispensing unit intended for formalin handling applications.

Design Characteristics of Controlled Dispensing Systems

Controlled dispensing stations typically incorporate:

- An enclosed dispensing area to limit vapour release during transfer
- Accommodation of bulk chemical containers within the enclosure
- Airflow management and vapour filtration to capture airborne contaminants
- A repeatable dispensing process to reduce operator variability



Role Within Risk Control Strategies

When assessed within the Hierarchy of Control, enclosed dispensing stations function as an engineering control, addressing exposure at source rather than relying solely on procedural compliance or PPE. They are commonly considered as part of a layered approach to chemical risk management.

Selection Considerations

Organisations considering controlled dispensing systems should evaluate:

- Compatibility with existing workflows
- Frequency and volume of dispensing activity
- Maintenance and filter management requirements
- Alignment with local COSHH assessments and safety objectives

Monmouth
Scientific

**PROTECTING
WHAT
MATTERS**

Controlled clean environments for critical applications.

Fume Containment | Biological Safety | Modular Cleanrooms | Powder Handling | Laminar Flow

UK Headquarters

Monmouth House,
Peninsula Business Park,
Bristol Road, Bridgwater,
Somerset,
TA6 4QB.

monmouthscientific.co.uk

info@monmouthscientific.co.uk
+44(0)1278 458090

