

# Circulaire® ATEX Fume Cupboards



**BRITISH  
MANUFACTURED**



**FIVE YEAR  
WARRANTY**



**EXPLOSIVE  
ENVIRONMENT  
RATED**



**ULTRA DEEP  
CARBON  
FILTERS**



The Circulaire® EX is an ATEX Zone 1 / T4 Rated Fume Cupboard range which incorporates the latest in containment technology to provide full compliance with COSHH regulations and BS7989:2001.

Compliance with the ATEX Essential Health and Safety Requirements is certified to EN13463-1: 2009.

The main filters consist of HEPA for particulate removal or Activated Carbon that will be specifically matched to your application and are monitored with a hydrocarbon detection system to assist in the replacement of the filters. Clean air is drawn through the front aperture by the fan at a speed of 0.50 m/s collecting any fumes that may be present before passing through the particulate pre-filter and carbon filter in the head section. The filtered air is then expelled back into the room.

## Size Specification (W x D x H)

EX 1100	EX 1400
1115mm x 823mm x 2176mm	1415mm x 823mm x 2176mm