

# Circulaire® Fume Cupboards



C 650

**BRITISH  
MANUFACTURED**



**FIVE YEAR  
WARRANTY**



**ULTRA DEEP  
CARBON  
FILTERS**



**CUSTOMISABLE  
FOR YOUR  
APPLICATION**



**The Circulaire® Recirculating Filtration Fume Cupboard is fitted with ultra-deep Activated Carbon or HEPA filters to guarantee maximum fume filtration efficiency.**

Fully compliant with COSHH regulations and international standards including BS7989:2001 for Filtration Fume/Particulate Cupboards.

To ensure the cabinet is functioning correctly, every model has an animated colour LCD display showing the airflow is working within safe and calibrated parameters.

All Monmouth products are designed to have both low noise levels and low energy consumption to help your company stay carbon neutral. LED Lighting further improve efficiency.

## Size Specification (W x D x H)

C 650	C 900
650mm x 650mm x 1125mm	900mm x 650mm x 1125mm